


# ***ENG 4550 – Introduction to Control Systems***

## ***Lab 6***



## **Lab 6: SRV02 Speed Control**

### **– Step response with PI control**

## 1. Lab report (Lab 6)

- Finish your lab report according to the template in Section 3.5.1 and tips in Section 3.5.3.

## II. RESULTS

Do not interpret or analyze the data in this section. Just provide the results.

1. Response plot from step 7 in Section 3.3.1.1, *Step response simulation with PI Control*
2. Response plot from step 9 in Section 3.3.1.2, *Step response implementation with PI Control*
3. Signal noise plot from step 10 in Section 3.3.1.2, *Step response implementation with PI Control*
4. Provide data collected in this laboratory (from Table 3.1).

# Submission of next lab

Section / Question	Description	Symbol	Value	Unit
Question 2	<b>Pre-Lab: PI Gains</b> Proportional Gain Integral Gain Open-Loop Time Constant Open-Loop Steady-state Gain	$k_p$ $k_i$ $\tau$ $K$		
Question 4	<b>Pre-Lab: DC Gain Estimate</b> DC Gain Estimate of $P_i(s)$	$ P_i(1) $		
Question 5	<b>Pre-Lab: Gain Crossover Frequency</b> Gain crossover frequency	$\omega_g$		
Section 3.3.1.1	<b>In-Lab: PI Step Response Simulation</b> Peak time Percent overshoot Steady-state error	$t_p$ PO $e_{ss}$		
Section 3.3.1.2	<b>In-Lab: PI Speed Control Implementation</b> Measured peak-to-peak ripple Steady-state error Peak time Percent overshoot	$e_{\omega, meas}$ $e_{ss}$ $t_p$ PO		
Section 3.3.2.1	<b>In-Lab: Step Response Simulation with Lead Control</b> Peak time Percent overshoot Steady-state error	$t_p$ PO $e_{ss}$		
Section 3.3.2.2	<b>In-Lab: Lead Speed Control Implementation</b> Peak time Percentage overshoot Steady-state error	$t_p$ PO $e_{ss}$		

- In 'ENG4550 control systems' on desktop, unzip 'Lab MatlabSimulink Software-20181001.zip' to a **NEW DIRECTORY**. All files you need in Lab 6 are in .../NEW DIRECTORY/Speed Control (Labs 6 and 7)
- When complete, **DELETE/REMOVE** your files and the **FOLDER** you created.



## 1. Configuring the SRV02 according to Section 3.4.2.

- In setup\_srv02\_exp03\_spd.m, make sure CONTROL\_TYPE is set to '**MANUAL**'. Run setup\_srv02\_exp03\_spd.m.

## 2. Follow the steps in **3.3.1.1 Simulation** .

## 1. Configuring the SRV02 according to Section 3.4.3.

- Setup q\_srv02\_pos.mdl: Double-click on the QUARC HIL Initialize block. Select the data acquisition device (q2\_usb or q8\_usb) you are using. Click on the **Defaults** and **OK** button.
- In setup\_srv02\_exp03\_spd.m, make sure CONTROL\_TYPE is set to 'MANUAL'. Run setup\_srv02\_exp03\_spd.m.

## 2. Follow the steps in 3.3.1.2 Implementing PI Speed Control.

- Before building the model (Step 7), click QUARC -> Set Default Options to avoid the possible target error.
- The amplitude in Figure 3.13 is wrong.